

Qualification of compressed gases
Microbiological Air Sampler for compressed gases
PDS Pharma III.Compressed.Gases.EX
Pre-programmed gases (compressed): Air, N2, CO2,
Ar and O2



### Description

The PDS Pharma III.Compressed.Gases.EX has been developed to control the microbiological contamination in compressed gases. The sampling method is based on the approved standard impaction principle but adapted to under pressure measuring scenarios. Standard 90-100 mm Petri dishes can be used. The most widely used gases in the industry are pre-programmed. All instruments are delivered as explosion-proof versions.

# Unique Benefit

The micro organisms are collected under the same pressure as in the gasdistribution system. The automatic decompressed cycle is started not before finish of the collection cycle. Pre-programmed also for Oxygen.

#### Application field

Our product has been designed for compressed gases used in the pharmaceutical and food industries and other critical fields. Compressed gases may be connected directly to the instrument from 1 up to 10 bars. The instrument is easy to handle, works absolutely automatically and is operated by a rechargeable battery. It is possible to connect it to a PC.

## Function al principle

The PDS Pharma III.Compressed.Gases.EX is a microbiological air sampler for compressed gases based on the principle of the Andersen air sampler ("golden standard") and aspirating air through a perforated plate. The resultant airflow is directed onto a standard agar poured plate. After the collection cycle, the Petri dish is incubated and the colony forming units (CFU) are counted. The system operates with a high-performance suction device, and the aspirated volume is continuously monitored. The system measures the inflow of air and regulates the aspirated volume to a constant value of 100 Standard litres per minute. The PDS Pharma III.Compressed.Gases.EX automatically regulates into solution according to pressure and temperature and thus gives all the time comparable results. For measurement is a functional call industry the recommended sampling volume — according to international guidelines — should be 100 litro.



## Accessories

Transport trolley (standard)
CG calibration unit (optional)
Additional sampling head (optional)

CAT Clean Air Technology GmbH Heimsheimer Strasse 35 D-70499 Stuttgart CEO: Steffen Hild Company site: Stuttgart Registered at destrict court Stuttgart (HRB18719) Bank details: Dresdner Bank in Stuttgart Account number: 1300 593 00 (BLZ: 600 800 00) Phone : +49 711 365 91 99 - 0 Fax : +49 711 365 91 99 - 99 Tax ID Number : 99 066 08216 VAT ID Number : DE 192 516 706